


# Applying Heat instructions



Osteopathy &  
Sports Injury Clinic

© Body Back-Up 2014



The hot or cold pack you heat in the microwave or cool in the freezer. The safe alternative for all ages to the hot-water bottle, gel packs and electric pads. No boiling water, trailing wires or danger from scalding.

The **Large / Medium Heat Pack** contains 100% naturally treated cereal encased in an inner cotton lining and comes with one outer slipcover.

Use your **Large / Medium Heat Pack** for comfort and relaxation while in bed, on the sofa, in the car or at work. The **Large / Medium Heat Pack** is ideal in helping to relieve:

**Tension and stress • Back pain • Stiff neck • Rheumatic pain • Menstrual pain  
Aching joints • Inflammation & bruising • Muscle strains  
Sciatic & migraine pains • Stomach aches**

The heat of the **Large / Medium Heat Pack** is slowly released and depending how and where you use it will stay warm for up to three hours snuggled up under blankets or a duvet, up to 45 minutes without insulation (perhaps on a stiff neck), up to an hour and a half (used as a back cushion). For the first few heatings you may experience variable heat retention and some condensation.

## • User Instructions •

**Heating:** The **Large / Medium Heat Pack** must **ONLY** be heated in a microwave on full power. Do not re heat the pack again until the pack has fully cooled from its use – should you require a second application of heat.

**650 – 750 Watts: up to three minutes • 800 Watts: up to two and a half minutes**

These times are for microwaves **with** a turntable. Microwaves **without** a turntable, your **Large / Medium Heat Pack** must be shaken and turned over halfway through the above heating times.

**Please note** we recommend the slipcover is removed before micro waving. To top up the heat pop it back in the microwave and heat for approximately **half the original time**.

**Freezing:** Remove the slipcover and place the **Large / Medium Heat Pack** in a plastic bag in the freezer for 2 – 4 hours depending on the setting of your freezer. Condensation may appear; gently shake and it will evaporate. When frozen the **Large / Medium Heat Pack** remains dry and supple and remains dry to the touch while defrosting.

## • Care Instructions •

The original channel design allows the contents to be positioned to the contours of your body. If the inner lining becomes dirty dab-clean with a damp wet cloth, do not immerse in water, dry slowly and thoroughly in a well ventilated place.

The outer slipcover is machine washable, see the instruction label. Remove all packaging before use and retain instructions for future reference.

Always follow the recommended heating times and **NEVER** over-heat. It is a natural product and not flame retardant.

[You may also be interested in our Fact Sheet — Applying Heat](#)

To find out more about Body Back-Up visit our website or email us at [practice@bodybackup.co.uk](mailto:practice@bodybackup.co.uk)

[www.bodybackup.co.uk](http://www.bodybackup.co.uk)

Info130 V4 (22.05.14)

For appointments call 01628 624544 or 020 8969 0247  
Maidenhead • London W9 (Maida Vale)